


| | |
|-------------------------------------|--|
| Project Name / Design | UK_IE_Calpol Vapour Plug & Nightlight_15ml_Leaflet |
| Date | 10.09.2020 |
| Country / Market(s) - Languages | English |
| Job / Version Number | 401995340_r3 |
| Technical Spec. / Drawing # | N/A mock up only |
| Component Type | Leaflet |
| Dimensions | 208mm x 185mm |
| New Component / Item Code | 86178 |
| Manufacturing / Packaging Fill Site | Zobebe, Italy |
| Old Component / Item Code | N/A |
| Printer Name | N/A |
| Pharma Code / Laetus Code | N/A |
| Print Method | N/A |
| Product / Eurocode / Mat. ID Code | 9216604 |
| Substrate | N/A |
| ArtWorks ID: EMEA_2020_00028511_000 | PC spec |
| | N/A |

ARTWORK SPECIFICS




REPRESENTATION: Colours represented with a diagonal line have been modified to aid PDF approval.

DISCLAIMER - It is the client's responsibility to analyze and review immediately upon receipt of digital files, films, plates or proofs prior to any recording of matrices and/or printing. In the event of a work has been found to be incorrect after the inspection described above, contact SCHAWK! to make sure that corrections are made before using this product. Schawk's liability is limited to fixing its own product. Schawk shall not be liable for any other damages and losses including but not limited to the costs of printing plates or cylinders, loss of print time, labor, substrates, paints, other materials and delays in delivery resulting from errors or defects in any of Schawk's products. Client is solely responsible for: (I) and failure on the file to comply with any applicable local, state or federal law or regulation, (II) and language errors in the digital file, proofs, plates, files or other media to properly advise the end consumer of any danger in the use or application of any product and, (III) and part of the file that violates the property rights of third parties acquired, without limitation, by the use, registration of copyright, patent or otherwise. Printer is responsible for apply and ensure barcode and datamatrix code readability

Kingsway North
Newcastle | Tyne & Wear
NE11 0JH | UK
+44 (0) 191 491 7777
schawk.com / www.sgkinc.com



| | | |
|------------|-------------------------|--------------------|
| Fonts used | HelveticaNeueLTStd-BdCn | Smallest Text Size |
| | HelveticaNeueLTStd-Cn | 7.0pt |
| | | Ave. Text Size |
| | | 12.0pt |
| | Largest Text Size | |
| | 12.0pt | |



Safety Information Sheet

Read this whole leaflet carefully before you use this product. Keep the leaflet, you might need it again.

Do not use with children under 3 months of age. The appliance can be plugged in and pads replaced by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.

Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision. The appliance is only to be used with the recommended vaporising medium. The use of other substances may give rise to toxic or fire risk. The appliance has hot surfaces to evaporate the active

ingredients and these surfaces should not be touched during use.

Do not touch the CALPOL Vapour Plug & Nightlight with wet hands or metal objects. Risk of electric shock. Use only in a wall mounted socket. Do not use with an adaptor or extension lead. Use with adequate ventilation. Do not use in small, confined areas. The vapour plug is not suitable for use in air conditioned environments. Do not place pads on varnished surfaces. Keep away from fabrics, bedding, and other materials that might present a fire hazard.

Revision date: August 2020.



Johnson & Johnson Ltd., High Wycombe, HP12 4EG, UK (Tel: 0808 238 9194)
Johnson & Johnson (Ireland) Ltd., Airton Road, Tallaght, Dublin 24, Ireland (Tel: 1800 22 00 44).
CALPOL® is a registered trade mark.

Directions For Use:

Cut across the top of the foil pouch to expose pad.

1. Using the foil to hold the pad, slide the pad into the slot of the vapour plug. When replacing a used pad, inserting a fresh pad will push out the old one. Dispose of used pads safely without touching them directly. Always wash hands after replacing pads.

2. Connect to power supply with the nightlight at the bottom. When using the vapour plug with a pad inserted, do not leave the unit in the plug socket for more than 16 hours at a time. After 16 hours, unplug the unit for a few hours before using again. Replace the pad before next use. Do not disassemble unit or nightlight. Nightlight goes on when the light goes out. When nightlight burns out replace entire unit. Do not block airflow of vapour plug unit or light sensor on bottom. If the light sensor is blocked the nightlight will stay on.

3. Unplug the vapour plug when not in use.

To Insert A New Pad

Allow the unit to cool for 10 minutes. Open a refill pad as described above. Inserting a new pad will push out the old one.

Be careful whilst removing and inserting the new pads as the device contains hot surfaces.

Remove and discard the old pad to use the nightlight alone.

Use only CALPOL® Vapour Plug Refill Pads with the CALPOL® Vapour Plug & Nightlight.